

YATLL INDUSTRY CO., LTD
YALU INDUSTRY CO., LTD



2024

PRODUCT CATALOGUE

Elderly care | Disabled care | Personal care



ABOUT US

YATTLL Industry Co., Ltd was founded in 2017, as one of the subsidiaries of YALU Industry Co., LTD, YATTLL is the proprietary brand of medical products.

YATTLL was co-founded by well-known Industry leaders and top Japanese engineers. It is committed to providing the society with more technological, convenient, and personalized products.

Through more innovative product design, smarter system solution, and more professional production experience, Yattll has always been committed to improving the quality of life and enhancing users' happiness.

YATTLL's product line includes two major directions: Smart travel and Rehabilitation assistance. The products cover multiple categories such as smart electric wheelchairs, smart scooters, smart hospital beds, and smart home-care beds. etc.

User's experience is always the first to YATTLL. YALU Industry Co., LTD set up its offices in Japan, Hong Kong and many cities in mainland China, which helps YATTLL not only do the deep researches of local market, but also collect feedback the first time. This makes YATTLL adjust its market strategy freely and offer better customer services in a very short time.

Our Mission:

To provide Innovative medical products and services with superior quality and competitive price, to improve the lives of consumers all over the world.



R&D: YATTLL focuses on Innovation, there are more than 10 experienced engineers and designers in our R&D team.

Marketing and Sales: We have set up branch offices in Japan, Hong Kong, Shenzhen, Huizhou. etc. Our products have exported to more than 30 Countries and areas.

Quality Control: Our company strictly complies with ISO quality management system. The quality control is stricter than ISO requirements, our implementation of production and testing standards have exceeded the industry standards.

DEVELOPMENT HISTORY

2003

Foundation of Yalu Industry Co., LTD

2015

Acquisition of a die casting company.

2019

Foundation of Micro switches company.

2021

Moving into YALU Techology Park.

2005

2005

Foundation of Home appliance company.

2016

Acquisition of an Injection molding company.

2017

Foundation of Yattll Industry Co., LTD

GLOBAL SALES

Hot -selling 126 countries and regions around the world



- | | | | | | |
|---------------|----------------|-----------|------------|---------------|---------------|
| 📍 UK | 📍 Saudi Arabia | 📍 Israel | 📍 Canada | 📍 Spain | 📍 Australia |
| 📍 Singapore | 📍 China | 📍 Ireland | 📍 Denmark | 📍 Finland | 📍 Netherlands |
| 📍 Dubai | 📍 Thailand | 📍 Germany | 📍 Malaysia | 📍 India | 📍 Brazil |
| 📍 New Zealand | 📍 Algeria | 📍 Bahrain | 📍 France | 📍 South Korea | 📍 Mexico |

OUR STRENGTH

Various large-scale production equipment and precision instruments



NEW



Model No.	YE246 II	YE200	YE200C	YE145D
Seat width	45.5 cm	49 / 52 / 56 cm	39 cm	45.5 cm
Seat depth	42 cm	43 cm	38.5 cm	50 cm
Front wheel	8" solid wheel	8"PU solid wheel	8"PU solid wheel	8" Suspension wheel
Rear wheel	12.5" solid, Suspension wheel	12.5"PU solid wheel	12.5"PU solid wheel	12.5" Solid, Suspension wheel
Battery	11Ah/15Ah/20Ah	6Ah/9Ah *2pcs	6Ah/9Ah *2pcs	11Ah/15Ah/20Ah
Motor	Brushless 250W	Brushless 250W	Brushless 200W	Brushed 250W
Charger	24V 2A/3A/5A	24V 2A/3A/5A	24V 2A/3A/5A	24V 2A/3A/5A
Controller	Brushless controller	Brushless controller	Brushless controller	Brushed controller
Brakes	Electromagnetic braking	Electromagnetic braking	Electromagnetic braking	Electromagnetic braking
Drive distance	15–18 km	15–18 km	15–18 km	10–13 km
Recline	NO	NO	Yes	NO
Weight capacity	135 kg	150 kg	135 kg	150 kg
Max speed	6 km/h	8 km/h	8 km/h	8 km/h
Net weight	25.5kg	28.8kg	32.8kg	26.5kg
Expanded Size	78*61*100 cm	98*60*110 cm	105*54*110 cm	107*61*100 cm
Page NO.	9P	10P	11P	12P

Best Seller

RECLINING BACKREST



Model No.	YE245C & MAX	Camel Hope YE246	YE245CR	YE200R
Seat width	45 cm	45 cm	45 / 52 / 56 cm	45 / 52 / 56 cm
Seat depth	43 cm	43 cm	43 cm	43 cm
Front wheel	8"Suspension wheel	8" Solid wheel	8"Suspension wheel	8" solid wheel
Rear wheel	12.5"Suspension wheel	12.5" Solid, Aluminum rim	12.5"Suspension wheel	12.5" Solid, Aluminum rim
Battery	11Ah/15Ah/20Ah	11Ah/15Ah/20Ah	11Ah/15Ah/20Ah	6Ah/9Ah *2pcs
Motor	Brushless 250W	Brushed 250W	Brushless 250W	Brushless 250W
Charger	24V 2A/3A/5A	24V 2A/3A/5A	24V 2A/3A/5A	24V 2A/3A/5A
Controller	Brushless controller	Brushed controller	Brushless controller	Brushless controller
Brakes	Electromagnetic braking	Electromagnetic braking	Electromagnetic braking	Electromagnetic braking
Drive distance	15-18 km	10-13 km	15-18 km	15-18 km
Recline	NO	NO	Yes	Yes
Weight capacity	135 kg / 150kg	135 kg	135 kg	150 kg
Max speed	8 km/h	6 km/h	6 km/h	6 km/h
Net weight	24.5kg	25.5kg	29.4kg	31.6kg
Expanded Size	98*62*104 cm	98*62*104 cm	98*62*103 cm	98*60*110 cm
Page NO.	13P	14P	15P	16P

HOT



Model No.	YE100	Camel Lite YE246	YEC35EBR	YE135	YE235
Seat width	49 / 52 / 56 cm	45.5 cm	40 cm	40/46/50/56/60 cm	40/46/50/56/60 cm
Seat depth	50 cm	43 cm	42.5 cm	42.5 cm	42.5 cm
Front wheel	8" Suspension wheel	7"solid wheel	8" Solid wheel	8" solid wheel	8" solid wheel
Rear wheel	12.5" Solid, Suspension wheel	8" solid wheel	12.5" Solid, Aluminum rim	24"steel wire wheel	12"steel wire wheel
Battery	11Ah/15Ah	11Ah/15Ah/20Ah	24V 12Ah	24V 24Ah	24V 24Ah
Motor	Brushless 250W	Brushless 250W	Brushed 250W	Brushed motor 200W	Brushed motor 24V 350W
Charger	24V 2A/3A/5A	24V 2A/3A/5A	24V 2A/3A/5A	24V 2A	24V 2A
Controller	Brushless controller	Brushless controller	Brushed controller	Brushed controller	Brushed controller
Brakes	Electromagnetic braking	Electromagnetic braking	Electromagnetic braking	Electromagnetic braking	Electromagnetic braking
Drive distance	15–18 km	15–18 km	10–13 km	18–20 km	18–20 km
Recline	NO	NO	NO	NO	NO
Weight capacity	150 kg	135 kg	135 kg	135 kg	135 kg
Max speed	6 km/h	6 km/h	6 km/h	8 km/h	8 km/h
Net weight	24.5 kg	21.6 kg	25.5 kg	42 kg	40.5 kg
Expanded Size	99*60*108 cm	98*58*99 cm	93*62*106 cm	94*65.5*113 cm	91*65*101 cm
Page NO.	17P	18P	19P	20P	21P

Meet the future



Model No.	Leader YE300	X RIDER	YECX1	Homecare Bed
Seat width	43 cm	43 cm	/	/
Seat depth	43 cm	43 cm	/	/
Front wheel	8"	8"	4"	5"
Rear wheel	8"	8"	4"	5"
Battery	9Ah * 2PCS	20Ah	5Ah	10Ah
Motor	Brushless 250W	24V 250W	/	/
Charger	/	/	/	/
Controller	Programmable Controller	Smart controller	Remote control	Remote control
Brakes	/	/	/	/
Drive distance	20 km	20 km	/	/
Recline	NO	NO	NO	Yes
Weight capacity	135 kg	125 kg	135 kg	150 kg
Max speed	6 km/h	6 km/h	Hand-propelled	Hand-propelled
Net weight	38 kg	20 kg	35 kg	130 kg
Expanded Size	88*55*120 cm	92*46*95 cm	80*55*110 cm	210*120*140 cm
Page NO.	22P	23P	24P	25P

YE246 - II

NEW

- A new generation of YE246.
- New designed backrest with extenable pocket on back.
- Fashion designed upholstery, removable and washable.
- New adjustable armrest, foldable for storage.
- Upgraded battery case, 1 click removable and load more.



YE200

Quick removable battery

- Intelligent programmable controller.
- Frame is made of high-precision craftsmanship.
- Back suspension system provides more outstanding experience.
- One step fold sequence.
- Quick removable battery design on two sides.
- Reinforced frame with a maximum loading of 150 kg.



12 Ah / 18Ah optional



YE200C

Special design for Children

- With multi-angle reclining backrest.
- Suspension system on the back frame.
- Waist-guard cushion offers better care to kids.
- Unique designed battery easily removable.
- Remote control for easier operating.



12 Ah / 18Ah optional



Anti puncture solid wheel



Suspension system on back frame



Soft & comfortable leg guard pad



Quick-released battery design

YE145D Brushed

- Net weight 28.4kg.
- 360 ° intelligent controller, easy to get started.
- Ergonomic-designed backrest, makes the seating more comfortable.
- Aircraft level aluminum alloy frame.
- It has dual high-performance motors 250W*2, strong power.

NEW



Suspension wheel,rubber material



Manual / electric mode switch



Aluminum alloy frame



One step fold

BEST
SELLER



YE245C & MAX

- Net weight 23kg, maximum load up to 135kg.
- Ergonomic designed cushion, removable & comfortable.
- Intelligent controller can be mounted on both the left and right sides.
- 11Ah / 15Ah / 20Ah Lithium battery available.

150 KG



CAMEL HOPE YE246

- It weighs 24.65KG(W/O) battery.
- Brushed motors 250W*2.
- Suspension on the back relieves the pressure.
- It has two unique dual base bar.
- Providing a bigger& wider seating space.
- With a small folding size.

BEST
SELLER



YE245CR

Reclining + Lightweight

- Lightweight electric wheelchair with Reclining backrest.
- With Medical motor and Cylinder to support smooth Reclining of backrest.
- Angle adjustable Headrest.
- Extendable and comfortable footrest.
- Loading capacity: 135KG



HOT
SELLING



Extended legrest



Intelligent universal controller



Front& rear perforated wheels



Brushless motor 250W each

YE200R

Comfort Max !

- Reclining + Wider Seating + 150KG
- Quick and easy removable battery, Lithium battery 6Ah X 2.
- Brushless motor 24V 250W each.
- Ergonomically designed cushion with reflects, super comfortable and removable.
- With Medical motor and Cylinder to support smooth Reclining of backrest.
- Angle adjustable Headrest.
- Extendable and comfortable footrest.
- Loading capacity: 150KG



Remote Control



Compact folding space



Electric/Manual Mode switch



Quick release battery design

HOT

YE100

Reinforced Frame

- Redesigned and reinforced frame which can load 150KG.
- Powerful & brushless motors 250W*2.
- With top advanced seat back cushion.
- It has a quick-reased battery.
- 8" rubber wheel, 12.5" rear wheel.
- Joystick can be mounted left and right side.



CAMEL LITE YE246

HOT

- Net weight 21.6kg,Max loading capacity: 135kg.
- With Top-quality brushless motors.
- Suspension on the front and back which can relieve stress easily.
- Folded in just ONE SECOND.
- It has a smallest turning radius 90cm.



Top-quality
brushless motor



Rear suspension system



Foldable footrest



Manual/ Electric
mode Switch



Backrest bag

YEC35EBR

Brushed

- Strong steel frame.
- Programmable controller, 360° free control.
- Flip up the armrest and foldable legrest.
- With an Electromagnetic braking system.
- With 12.5" rear wheel, Aluminium rim.
- With Lead Acid batteries 24V 12Ah.



HOT



YE135

Big Rear Wheel



- Programmable controller and 360° precise control.
- Comfortable and soft cushion, detachable for cleaning.
- Controller can be mounted on Left or Right of the chair.
- 24" rear wheel makes the chair drive faster.
- More wheel options available.



YE235

Portable Wheelchair

- Classic design with small size.
- Optional seat width.
- More controller and battery options.
- The upgraded armrest has adjustable height and length.
- Simple folding.



Armrest move forward and backward



Elevating armrest



Foldable backrest



Swing away footrest

Leader YE300

Comfort Scooter

- One press auto folding.
- McNamee Wheel makes the turning freely.
- Easy removable battery, 12Ah / 18Ah optional.
- Comes with comfortable headrest, height adjustable.
- Storage basket on back of wheelchair.



X RIDER

Smart Lightweight Scooter

- LED screen with multi function display.
- Magnetic charging house for cell phone.
- SOS alarm with location tracking system.
- Compact design with a small folding size.
- Automatic headlights.



YECK1

Transfer Chair

- Lightweight and foldable transfer chair.
- Armpit support to simulate a pick-up motion.
- Aluminum frame with a unique folding way.
- Height adjustable. Remote control.




Homecare Bed

- Smart body condition monitoring
- Multi-angle turnover
- Electronic display screen
- Enhancement size of the care bed
- making it suitable for home use



Electronic display screen



Migrant wheel



Intelligent monitoring



Manual Wheelchair



YE119 Steel wheelchair

- Steel wheelchair, chrome plated.
- Fixed arm and fixed leg.
- 6" Caster with PVC tyre.
- 24" Spoke wheelchair with solid tyre.
- Double cross bar.
- N.W: 13.3kg.
- Size: 104*68*87.5cm.



YML120 Steel Wheelchair

- Strong steel wheelchair.
- Fixed arm and fixed leg.
- Lightweight and easy foldable.
- 7" Caster with PVC tyre.
- 16" Plastic rear wheel.
- N.W: 12kg.
- Size: 80*65*90cm.



YML121 Lightweight aluminum wheelchair

- Material: 2.0 mm aluminum.
- Flip up armrest with PU pad.
- Reliable braking system.
- Foldable footrest.
- 6" caster, 8" solid rear wheel.
- Seat height: 42cm.
- N.W: 7.5kg.

Manual Wheelchair



YE122 ULTRA Lightweight wheelchair for Children

- Aluminum frame.
- Flight free wheelchair.
- 6" castor, 8" rear wheel.
- Seat width 34 cm, seat height 50 cm.
- Fold and go.
- N.W: 7.5KG.



YE123 Steel wheelchair for Children

- Fixed armrest and fixed legrest.
- 7" castor, 20" mag wheel.
- Durable mesh cushion, detachable, washable.
- Handle brake and rear wheel brake.
- N.W: 11KG.



YE126 Aluminum frame

- Flip up armrest, foldable legrest.
- Fold down backrest to save storage spaces.
- 6" castor, 20" mag wheel.
- Durable mesh cushion, detachable, washable.
- Handle brake and rear wheel brake.
- N.W: 10KG.

Wheelchair parts

Motor Pair / Remote control



YSPM005
320W D.C Motor



YSPM001
250W Brushed Motor



YSPM004
250W Brushless Motor



YSPM002
250 Brushless Motor



YSPCR001
Remote control

Controller



YSPC002
Controller (brushless)



YSPC003
Controller (brushless)



YSPC001
Controller (brushed)



YSPC004
Controller (Reclining)



YSPC005
Controller (LED screen)

Battery




YSPBL360
6Ah/9Ah Quick removable



YSPBL210
11Ah/15Ah Quick removable



YSPBL110 / 150 / 200
Lithium battery



YSPBL401
26Ah One click removable



YSPLA11
12Ah Lead acid battery

Wheelchair parts

Wheels / Charger



YSPW1202
12.5" PU Solid wheel



YSPW1204
12" Suspension wheel



YSPW1205
12.5" PU Solid wheel



YSPW1206
12.5" PU Solid wheel



YSPW701
7" PU Solid wheel



YSPW801
8" PU Solid wheel



YSPW803
8" PU Solid wheel



YSPW804
8" Heavy duty wheel



YSPCG212
YSPCG313 / YSPCG515
2A/3A/5A Charger



YSPCG203
2A Charger for
Lead acid battery

Accessories



YSPSU001
Foamed upholstery



YSPSU002
Foamed upholstery



YSPSU002H
Heating upholstery



YSPSU001H
Heating upholstery



YSPBG12
Storage bag

Wheelchair parts

Accessories



YSPBG10
Travel bag



YSPBG11
Armrest sidebag



YSPAR001
Armrest assembly (L&R)



YSPFR001
Footrest assembly



YSPFR002
Footrest assembly



YSPAC001
Bracket for Rear Controller



YSPAC002
Controller holder (short)



YSPAC002
Controller holder (long)



YSPWK52
YSPWK56
Seat width widening kit



YSPCH001
Cup holder



YSPAR001-A01
PU arm pad



YSPAR001-A02
Curved PU arm pad



YSPPH002
Cup & Phone holder



YSPPH001
Phone holder



YSPHR001
Headrest, Angel adjustable



YATLL INDUSTRY CO., LTD
YALU INDUSTRY CO., LTD

Tel: +86 755 81707199

Email: info@yattll.com

Website: www.yattll.com

Factory: Yalu Technology Park, No. 88 Pingshan Industrial Road,
Tangxia Town, Dongguan, 523728, Guangdong, China

Office: Rm 907, A&I International Centre, No.8 Gaoerfu blvd,
Longhua Dist, Shenzhen, Guangdong, China

